

Analytik Jena UK Open For Business At L.A.B Exhibition, ExCel Centre, London, 2-4 October 2007

We are delighted to announce that **Analytik Jena**, a world leader in the manufacture of high quality analytical instrumentation for a wide variety markets, **has opened a UK office**. The L.A.B Exhibition at the ExCeL Centre, London 2-4 October, will provide the first opportunity to meet the AJ-UK team and to see some of the superb equipment. **Come and see us at Booth D09.**

On show at the exhibition will be the multi EA[®] 3100 and the novAA[®] 300.

The **multi EA[®] 3100** is a new generation, top of the range system for fully automatic **simultaneous determination of Total Nitrogen and Total Sulphur** with sequential determination of Total Chlorine in solid, liquid, paste, viscous and gaseous samples. **The instruments are especially well suited for the analysis of petrochemical and related samples, such as oil, fuel, solvents, polymers, gas and liquid gas analysis**, as well as fine chemical, EOX and AOX analytics. With automatically switching detection methods for each element, detection limits of 10 ppb S, 50 ppb N and 0.1 ppm Cl can be achieved.

The **novAA[®] 300** is a fully automated entry-level flame atomic absorption **spectrophotometer with features and performance normally found only on much more expensive systems**. High quality optics, computer controlled gas box, a motorized lamp turret and automatic control of flame conditions ensure high quality results are obtained quickly and easily. With an automatic 6-lamp hollow cathode lamp turret and D2 HCL background correction, accurate analyses can be carried out even on samples with high background content.

We hope to see you at the show, but if you can't make it to the show, don't worry – **our brand new website, www.analytikjena.co.uk, has details of the multi EA[®] 3100 and the novAA[®] 300 and the rest of our product range**. If you'd like any **more information on AJ-UK**, please contact Vince Phelan or David James on 024 7647 0034 or sales@aj-uk.co.uk.